

# National Laboratory Service

QA/QC

David Johns

**create a better place**

# National Laboratory Service

## Summary

Method is UKAS accredited to ISO17025

Analyse CRM (BCR No.130) with samples

Mastersizer serviced annually

Sieves are certified

Balances calibrated annually and checked daily

Three measurements taken for each sample

**create a better place**

# National Laboratory Service

## Advantages

Provides confidence in reported data

**create a better place**

# National Laboratory Service

## Disadvantages

Difficult to obtain CRM that covers full range of band-sizes

**create a better place**

# National Laboratory Service

The End

**create a better place**